


70+ Interview Questions for Pharmaceutical Microbiology

/ Pharma Articles / By tech-publish.com

Top 70+ Interview Questions for Pharmaceutical Microbiology



This page covers most of the interview questions and answers asked during a technical interview round of Pharmaceutical Microbiology.


You will find interview questions and answers on microbiology basics, pharmacopoeial chapter and sections applicable for microbiology, Culture Media, Growth Promotion test, microbiological quality control for sterile and non-sterile dosage forms, good practice in a microbiological laboratory, aseptic technique to be followed in the microbiological laboratory, pure culture, Pyrogens and endotoxins, different staining techniques used in the pharmaceutical microbiology laboratory, microbial identification, Gram-staining technique, fungal staining, approved culture collections, sterility testing and environmental monitoring.

The interview questions cover questions from basic to advance level of technical aspects. These in

X



1. Explain the parts of compound light microscopes used in the pharmaceutical microbiology laboratory?



- **Illuminator:** Light source at the base of the microscope enable better view of object;
- **Condensor:** This is a two lens system that collects and concentrates light from the illuminator and directs it to the iris diaphragm;



- **Iris diaphragm:** regulates the amount of light entering the lens system;
- **Stage:** This is a platform for the slide with hole in the center to let light from the illuminator pass through. It consists of clips to hold the slide lace;
- **Eye pieces:** To view object;
- **Objectives:** To magnify the object as per rating on it

2. What are the pharmacopoeial chapter and sections applicable for microbiology?

Pharmacopoeial chapter(s)	Use
Microbiological examination of nonsterile products: total viable aerobic count (Ph. Eur. 2.6.12, USP <61>) USP <2021> Microbial enumeration tests-nutritional and dietary supplements USP <2023> Microbiological attributes of nonsterile nutritional and dietary supplements	Organism count in raw materials, water, finished products
Microbiological examination of nonsterile products: tests for specified organisms (Ph. Eur. 2.6.13, USP <62>) USP <2022> Microbiological procedures for absence of specified microorganismsnutritional and dietary supplements	Type of organisms present in raw materials, water, finished products
Ph. Eur. 5.1.4 Microbiological quality of pharmaceutical preparations/USP <1111> Microbiological attributes of nonsterile pharmaceutical products	Determining limits and control factors



Antimicrobial preservative efficacy testing (Ph. Eur. 5.1.3, USP <51>)	Antimicrobial preservative efficacy
Microbiological assay of antibiotics (E.P 2.7.2., USP<81>)	Assays of antibiotic
USP <1112> Application of water activity determination to nonsterile pharmaceutical products	Water activity can be used to predict microbial proliferation in the product
USP <1211> Sterilization and sterility assurance of compendial articles. There are a series of subchapters that describe specific sterilization methods	Sterilization and microbial reduction
USP <55> Biological indicators	Biological indicators for assessing microbial reduction
USP <1113> Microbial characterization, identification, and strain typing	Microbial identification
USP <1117> Microbiological best laboratory practice	Management of microbiological good laboratory practice
USP <1116> Microbiological control and monitoring of aseptic processing environments	Environmental monitoring for aseptic environments
Ph. Eur. 5.1.6 Alternative methods for control of microbiological quality/ USP <1223> Validation of alternative microbiological methods	Alternative methods/ Rapid microbiological methods
USP <1115> Bioburden control of nonsterile drug products	Bioburden control
USP <1227> Validation of microbial recovery from pharmacopeial articles	Microbial method validation

3. Which factors of culture media affects the cultivation of microorganisms?

Culture media factors affecting the cultivation of microorganisms are optimum nutrients, oxygen or other gases, moisture, pH, and temperature.



4. What are the important nutrients of culture media that affects the cultivation of microorganisms?

The key nutrients of culture media that affects the cultivation of microorganisms are sources of carbon, nitrogen, water, inorganic phosphates and sulfur, vitamins, and trace metals.

5. Whose name recognize as “The Father of Culture Media”?

Robert Koch (1843–1910) discovered that broths based on fresh beef serum or meat extracts (so-called bouillons, the term “broth” for liquid culture medium being analogous to broth or soup) produced optimal growth.

6. Who was the first scientist cultivated the microorganisms on a growth medium?

The French chemist and microbiologist, Louis Pasteur (1822–1895).

7. Petri dish is named based on which scientists' contribution?

Types of medium	Use
Nutrient agar or broth and tryptone soya agar or broth. Tryptone soya agar (equivalent to soyabean casein digest medium)	→ Environmental monitoring→ Isolation and cultivation of nonfastidious and fastidious microorganisms Note: Fastidious bacteria are bacteria that need special nutritional supplements and conditions to grow on agar plates. Nonfastidious bacteria are bacteria that do not need special nutritional supplements and conditions to grow on agar plates
Tryptone soya broth	Sterility testing and as a general growth broth in microbial enumeration tests, as well as used for media simulation trials
vegetable peptone broth	Media filling trials
fluid thioglycollate medium	used for the growth of bacteria (aerobic and anaerobic) as a part of the sterility test
fungi	Sabouraud dextrose agar or malt extra agar
R2A	Microbiological examination of water (This is a low nutrient agar used for the cultivation of heterotrophic microorganisms).
Columbia blood agar	Detection of hemolytic reactions by Staphylococci

9. What is the process of Growth Promotion test of microbiology culture media? What should be the acceptance criteria?

1. Medium to be inoculated with a microorganisms <100 Colony Forming Units (CFU).
2. Use not less than five unique strains as recommended by the pharmacopoeias and two organisms from environment isolates on rotation basis (It should be within the five passages from the original reference culture seed lot).
3. Compare the growth on the medium with medium previously approved culture lot.
4. Acceptance is no more than a factor of 2 differences in productivity ratio calculations.

10. What is the calculation of productivity ratio and what is the acceptance criterion for agar media plate or solid media for growth promotion test?

Productivity ratio =



Quantitative techniques of Growth Promotion Tests are of two types.

- Ecometric
- Miles – Misra (Drop Count)

12. Explain ecometric method of Growth Promotion Tests?

This technique is semi-quantitative variant of the streaking.

This technique is operator-dependent and has lower precision. However, it can be used to great effect with practice.

Following are the steps:

- A fresh suspension of the challenge organism is taken into a calibrated loop (One loopful of inoculum).
- Five streaks are streaked out into four quadrants onto the agar plate along with a final streak in the center of the plate.
- These plates are then incubated overnight for growth.
- The patterns of growth are interpreted to provide an Absolute Growth Index (AGI)

13. Explain Miles – Misra technique (the drop count technique) of Growth Promotion Tests?

This technique involves spreading droplets of known quantities (10 µL) of microbial suspensions.

The test plate is compared against control plate after incubation to verify number of colonies recovered.

Note:

The accuracy of the method is dependent on following factors:

- Dilution used
- Number of colony forming units (cfu) in the inoculum
- Volume of the inoculums
- Spreading technique



- Copious Growth
- End-point Methods
- Most Probable Number (MPN)
- Kinetic Parameters

15. Explain Copious Growth method of Growth Promotion Tests for Broth Media?

This is the method of choice as per compendial. The challenge of broth media is done by comparing the growth over the period of time with a control batch (which provides a qualitative assessment of copious growth). In this method, the challenge organism is inoculated at a very low level (< 100 CFU per unit) and incubated at the prescribed temperature for the prescribed period of time (3 days or 5 days).

The advantage of this method is that it does not require a great deal of labor. Semi-quantitative assessment can be done by constructing a growth index from slight to copious growth (normally a scale of +, ++, or +++).

16. Explain End-point Methods method of Growth Promotion Tests for Broth Media?

In this method, very low levels of inoculum are added to multiple tubes of the two media under testing. Growth frequency is compared between the two media to understand the equivalency.

17. Explain Most Probable Number (MPN) method of Growth Promotion Tests for Broth Media?

This is a Microbial Limits Test. In this method, the unknown sample is prepared in a ten-fold dilution series and added to nutrient broth in replicate tubes. The tubes will then either turn turbid (growth) or remain clear, and allow for an estimate of the most probable number of microorganisms.

18. Explain Kinetic Parameters method of Growth Promotion Tests for Broth Media?



19. How growth promotion tests for selective media is performed.

Selective media requires a different approach than general purpose media. The objective of the test for selective media is to see that the media supports the growth of specific microorganisms. Inhibition, colony morphology and pigment are verified during the testing. In order to determine growth of selective media, specific microorganisms are used as positive or negative indicators.

20. What are the steps involved in the manufacturing of culture media?

a. Collection of prerequisite

- Decide batch size
- Vessel for holding or dispensing the media
- Weighing balance
- Petri dishes/ Containers
- Dehydrated media

b. Rehydration of media

- Media powder is rehydrated by mixing the medium in the volume of water (Amount should be based on manufacturer's instruction)
- Homogenize the solution by mixing
- Care must be taken to avoid scorching the media.
- The media should clarify near boiling (95–100 °C), and the media should only be allowed to boil less than 1 min.

c. Sterilization

d. Addition of supplements (depending on the types of media requirement)

e. Filling in container

f. Status labeling

21. What are the acceptance criteria for microbiological quality of non-sterile dosage forms and raw materials?



Aqueous preparations for oral use	100	10	Escherichia coli absent in 1 g or 1 mL
Rectal use	1000	100	Need based
Oromucosal use Gingival use Cutaneous use Nasal use Auricular use	100	10	Staphylococcus aureus absent in 1 g or 1 mL Pseudomonas aeruginosa absent in 1 g or 1 mL
Vaginal use	100	10	Staphylococcus aureus absent in 1 g or 1 mL Pseudomonas aeruginosa absent in 1 g or 1 mL Candida albicans absent in 1 g or 1 mL
Inhalation use	100	10	Staphylococcus aureus absent in 1 g or 1 mL Pseudomonas aeruginosa absent in 1 g or 1 mL Bile-tolerant gram-negative bacteria absent in 1 g or 1 mL
Transdermal patches	100	10	Staphylococcus aureus absent in 1 g or 1 mL Pseudomonas aeruginosa absent in 1 g or 1 mL
Raw material for pharmaceutical use	1000	100	Need based

22. What are the risks because of microorganism in the pharmaceutical formulation if not within the acceptable limit?

- Decomposition of the product
- Degradation of product
- Infection to the patient

23. Give few examples of good practice in a microbiological laboratory.

- Required equipment and media should be sterilized prior to use.
- Sterilized equipment and media should be maintained and stored in such a way that it should not be contaminated.
- Frequently disinfect hands and working surfaces while carrying out aseptic operations.
- Fly, rodent and pest control.



- Sterilize contaminated used media and waste prior to disposal.
- Pipetting should not be done mouth.

24. Why aseptic technique should be followed in the microbiological laboratory?

To prevent contamination and cross-contamination in microbiology laboratory. Probability of contamination in the laboratory is – samples, media, environment, facility etc.

25. How asepsis can be achieved in the laboratory?

- Frequent washing and disinfecting hands.
- Use of unidirectional airflow station for critical operations.
- Handling of positive control at the end of analysis.
- Not touching the samples and accessories directly by hand.

26. What is pure culture?

A pure culture is a population of cells or multicellular organisms growing in the absence of other species or types in which cells are genetic clones of one another.



27. What is most common method to isolate individual cells and produce a pure culture?

The most common method to isolate individual cells and produce a pure culture is to prepare a streak plate method.

28. Explain Streak Plate Method.

- a. Sterilize the inoculating loop by heating it until red hot in a flame until it is red hot. Allow it to cool.
- b. Pick up a loop full of liquid inoculum or bacterial growth from the surface of an agar plate and, starting about 2.5 cm in from the edge of the plate, streak lightly back and forth with the loop flat, making close, parallel streaks back to the edge of the plate to a quarter of the plate.
- c. Sterilize the loop and cool again, going back to the edge of area first quarter that you just streaked, extend the streaks into the second quarter of the plate.



d. Repeat the same process for remaining two quarter of the plate.

e. Flame and cool the loop again.

29. What should be done during microbiology testing if sample exhibit antimicrobial activity?

When the sample possesses antimicrobial activity, it that requires neutralization. Following are commonly used principle for neutralization:



30. What are pyrogens and endotoxins?

Pyrogens and endotoxins are a heterogeneous group of chemical entities that cause fever when injected.

31. What is the source of pyrogens and endotoxins?

Pyrogens can be nonbacterial as well as bacterial in origin.

In the pharmaceutical industry, mostly observed source is from Gram-negative bacteria. That is the lipopolysaccharide (LPS) from the bacterial cell wall.

Pyrogens are the metabolic byproduct of the Gram-negative bacteria.

32. Explain types of microbial identification methods.

Identification methods can be divided into two types:

- Phenotypic
- Genotypic

33. What are the differences between Phenotype and Genotype?

- “Genotype” is an organism’s full hereditary information.
- “Phenotype” is an organism’s actual observed properties, such as morphology, development or behavior.

34. What are the different staining techniques used in the pharmaceutical microbiology laboratory?

- Gram-stain
- Bacterial spore stain
- Fungal staining
- Ziehl-Neelsen stain

35. What is the objective of microbial identification?

The objective of microbial identification is to differentiate one microbial isolate from another with re



• Genus: group of similar species

- Species: group of similar strains
- Type: strains within a species
- Strain: an isolate of a particular species

37. What is the first step of microbial identification?

The first step of most identification program is colony and cellular morphology of the microorganism.

Colony morphology is a method that used to describe the characteristics of an individual colony of bacteria on agar media in a Petri dish.

Colony morphology is normally classified based on the form, elevation and margin. This can be further classified as follows.

X



39. How the Gram-staining technique is useful?

Gram-staining technique that allows to visualize the morphological types of bacteria using a compound light microscope under magnification of 100x.

40. Explain the principle of the Gram-staining technique?

Step 1: Crystal violet (primary stain) aqueous solutions consist of +ve and -ve ions. These ions penetrate through the cell wall and cell membrane of all types of bacteria. The +ve ion stains the bacterial cells and stains the cells purple.

Step 2: Iodide (the mordant) interacts with crystal violet and forms complexes.

Step 3: Decolorizer (made of acetone and alcohol 95%) interacts with the lipids of the cell membrane. A Gram-negative cell loses its outer lipopolysaccharide membrane, and the inner peptidoglycan layer is left exposed.

On the other hand, a Gram-positive cell becomes dehydrated because of property of decolorizer. The complexes trapped within the Gram-positive cell due to the multilayered nature of the peptidoglycan.

Step 4: When decolorization is added, the Gram-positive cell remains purple, and the Gram-negative cell loses its purple color.

Step 5: A Safranin, the counterstain is used to provide color to Gram-negative bacteria a pink/ red color.



41. What is the procedure of Gram-staining technique?

- a. Take a clean, grease free slide.
- b. Prepare the smear of suspension on the clean slide with a loopful of sample.
- c. Air-dry or heat-fix smear of cells for around 1 minute with crystal violet staining reagent. (**Note:** Too heavy or too light cell concentration will affect the Gram Stain results.)
- d. Rinse the slide with gentle stream of tap water for 2-3 seconds.
- e. Flood the gram's iodine for 1 minute and wash with gentle stream using tap water for 2-3 seconds.
- f. Flood slide with 95% alcohol or acetone for about 15-20 seconds until decolorizing agent running fr



i. Blot dry with absorbent paper.

j. Observe under oil immersion using a Brightfield microscope.

k. Results:

Gram-negative bacteria will stain pink or red

Gram-positive bacteria will stain blue or purple

42. What is the procedure of Gram-staining technique?

a. Primary Stain: Crystal Violet Staining Reagent.

Solution A for crystal violet staining reagent

Crystal violet (90% dye content) – 1g

2g Ethanol, 95% (v/v) – 10 ml

Solution B for crystal violet staining reagent

Ammonium oxalate, 0.4 g

Distilled water, 40 ml

X



Mix solution A and solution B to obtain crystal violet staining reagent. Store for 1 day and then filter using filter paper.

b. Mordant: Gram's Iodine

Iodine, 0.5 g

Potassium iodide, 1.0 g

Distilled water, 150 ml

Triturate iodine and potassium iodide in a mortar and add water slowly with continuous trituration until the iodine is dissolved. Store in amber bottles.

c. Decolorizing Agent

Ethanol, 95% (vol/vol)

Alternate decolorizing agent – 1:1 acetone and ethanol mixture.

d. Counterstain: Safranin



Use 1.25g Safranin O and mix it with 50 ml 95% Ethanol (Solution A). Take this Solution A 5 ml and mix it with 45 ml distilled water

43. Which indicators are used for spore staining?

Malachite green (a triarylmethane dye) and a safranin (an azonium compound) counterstain is useful tool in identifying the presence or absence of spores.

This is referred to as the Schaeffer-Fulton stain.

44. Explain the Schaeffer-Fulton staining method for endospore staining technique?

1. Prepare a bacterial smear on a clean slide, air dry and gently heat fix.
2. Cover the slide with a piece of paper towel, and place on a staining rack over the water bath with boiling water.
3. Flood the paper towel on the slide with Malachite green (primary stain). Steam the slide for about 5 minutes.
4. Remove the slide from the water bath, and remove the paper towel from the slide.




5. Examine the slide with a light microscope under an immersion objective

Result:

Endospores appears green

Vegetative cells appear red or pink



45. What is Lactophenol cotton blue stain or fungal staining?

Lactophenol cotton blue stain is used to examining yeast and filamentous fungi under microscope.

Phenol acts as a disinfectant to kill living organisms


Lactic acid used to preserve the fungal structures

Cotton blue is used to stain fungal cell wall and other fungal structures

On staining, it fungus will be with blue colored spores and structures, such as hyphae.

The identification of fungi using macroscopic and microscopic techniques is difficult and requires a tr.





46. What reagents used to prepare Lactophenol cotton blue staining solution?

Distilled water

Cotton Blue or Aniline Blue

Phenol Crystals

Glycerol

Lactic acid

70% ethanol

47. What is use of Ziehl-Neelsen stain?

The objective of Ziehl-Neelsen stain is to differentiate bacteria into acid fast group and non-acid fast groups.


48. What Acid-Fast stain called as Ziehl-Neelsen stain?

This technique was first developed by Ziehl and later on modified by Neelsen, therefore, this method is also called as Ziehl-Neelsen staining techniques.





Reagents	Color of Acid fast	Color of Non-acid fast
Primary dye – Carbol fuchsin	Red	Red
Decolorizer – Acid alcohol	Red	Colorless
Counter stain – Methylene blue	Red	Blue



49. What is Total viable aerobic count?

Total viable aerobic count is designed to count the number of microorganisms (as colony forming units, CFUs) in a non-sterile pharmaceutical product or raw materials.



is further divided into two parts, total aerobic microbial count (TAMC), and total yeast and mold counts (TYMCs).

50. Which are the methods used for analysis of total viable aerobic count?

- Membrane filtration technique – Sample is filtered and the filter is placed on defined media
- Pour plate technique – Sample aliquot is taken and placed in a Petri dish and specified media poured onto the sample
- Spread plate technique – Sample aliquot is placed on the surface of defined media and smeared evenly over the surface
- Most probable number (MPN) technique – This method is used for mainly insoluble materials. The sample dilutions are placed into a series of replicate tubes and the number of tubes showing growth give a statistical evaluation of the number of microorganisms in the sample.

51. Which is Bioburden test Method Validation?

The Bioburden Method Validation is done to demonstrate the adequacy of sample preparation method and the ability of the media to recover microorganisms in the presence of the test sample. Following should be considered during the method validation.

(a) Growth promotion of media

X

(b



- American Type Culture Collection (ATCC)
- National Collection of Industrial and Marine Bacteria (NCIMB)
- Collection of Institute Pasteur (CIP)
- Imperial Mycological Institute (IMI)
- National Collection of Pathogenic Fungi (NCPF)
- National Biologicals Resources Centre (NBRC)



53. Explain the sample preparation process for Bioburden test.

- **Sample handing area:** Sample should be prepared on a laboratory bench, within unidirectional airflow cabinet or an isolator.
- **Bacteria cultures handling area:** Biosafety cabinet or Microbiological Safety Cabinet should be used.
- **Sample preparation for Water-soluble products:**

a. Dissolve or dilute the product with ration of 1 in 10 dilution in phosphate buffer solution pH 7.2, If necessary, adjust to a pH of 6 – 8.

b.



- **Sample preparation for Non-fatty and insoluble material in water:**

- a. Suspend the product with 1 in 10 dilution in phosphate buffer solution pH 7.2.
- b. A surfactant such as 1 g/L of polysorbate 80 can be used to assist the suspension of poorly wettable substances.
- c. If require, adjust to a pH of 6–8.
- d. If required, further dilutions can be done to get not more than 250 CFU/plate in case of TAMC, 50 CFU/plate in case of TYMC.

- **Sample preparation for fatty products:**

- a Dissolve in isopropyl myristate sterilized by filtration, or mix the product to be examined with the minimum necessary quantity of sterile polysorbate 80.
- b. Heat if required for NMT 40 °C or, in exceptional cases, to not more than 45 °C and maintain the temperature in a water bath.

- c. Add pre-warmed diluent to make a 1 in 10 dilution of the product.
- d. Form of an emulsion.

X

c.



A spectrophotometer method is used to measure turbidity of medium suspension. When light passes through a suspension of microorganisms, light gets scattered and the amount of scatter is an indication of the concentration present in the suspension.

While doing the estimation using this method, a calibration curve needs to be constructed using series of the known concentration. Based on the calibration curve, unknown concentration can be easily identified.

55. What are the common sources of pyrogens?

Non-bacterial source:

- Antigens
- Poly nucleotides
- Steroids
- Adjuvants
- Viruses
- Fungi

Bacterial source:

- Streptococcal toxins
- Staphylococcal enterotoxins



- Bacterial cell wall – lipopolysaccharides (endotoxins)

56. What are the method used to test pyrogens and endotoxins?

- Rabbit pyrogen test
- Limulus Amoebocyte Lysate (LAL) testing for bacterial endotoxin

57. Is sterility test qualitative or quantitative test?

The sterility test is a qualitative test. Results are written as presence or absence of turbidity based on growth for bacteria and fungi in media.

58. Which are commonly used sterility testing methods?

- Membrane filtration: This is further classified as open method and closed system method.
- Direct inoculation method

59. What components required performing sterility testing by membrane filtration method?

- Filter with pore size 0.45 μm
- Filter diameter about 50 mm
- Cellulose nitrate filters for aqueous, oily, and weak alcoholic solutions
- Cellulose acetate filters for strong alcoholic solutions

60. What are the growth mediums used to perform sterility testing?

- Fluid Thioglycollate Medium for anaerobic bacteria and it also isolates aerobic bacteria
- Soya-bean Casein Digest Medium for the isolation of fungi and aerobic bacteria

61. What types of growth promotion test (GPT) to be done for growth medium used to perform sterility testing?



- GPT for Fluid thioglycollate medium shall be done using *Clostridium sporogenes*, *P. aeruginosa*, *S. aureus*;
- GPT for Soya-bean casein digest medium shall be done using *A. niger*, *B. subtilis*, *C. albicans*.

The medium shall be inoculated with not more than 100 CFU and incubated for 3 days (bacteria) or 5 days (fungi).

Result: Clearly visible growth must be observed.

62 What should be incubation conditions and times for sterility testing?

Incubation conditions and times for sterility testing are as follows:

- Fluid Thioglycollate Medium at 30 – 35 °C
- Soya-bean Casein Digest Medium at 20 – 25 °C
- Total incubation period of 14 days with visual examination for turbidity.

63. How to interpret sterility test validation outcome?

There are two outcomes with the sterility test:

- (1) Clearly visible growth – test sample and the control tubes are equivalent
- (2) No clearly visible comparable growth – Test needs to be modified and the validation needs to be repeated

64. What should be the environment for sterility testing?

The sterility test environment should be EU GMP Grade A with a Grade B background or a Grade A isolator operator with background Grade C or D.

X

64



to be collected when unplanned intervention occurred during the batch.

66. What is the composition of Fluid thioglycollate medium?

The composition of Fluid thioglycollate medium is as follows:

Name of ingredient	Weight
L-Cystine	0.5 g
Agar	0.75 g
Sodium chloride	2.5 g
Glucose monohydrate/anhydrous	5.5/5.0 g
Yeast extract (water-soluble)	5.0 g
Pancreatic digest of casein	15.0 g
Sodium thioglycollate or	0.5 g
Thioglycollic acid	0.3 ml
Resazurin sodium solution (1 g/l of resazurin sodium), freshly prepared	1.0 ml
Water R	1000 ml

pH after sterilization 6.9 to 7.3.

67. What is the composition of Soya-bean casein digest medium?

The composition of Soya-bean casein digest medium is as follows:

X



Papain digest of soya-bean meal	5.0 g
Sodium chloride	5.0 g
Dipotassium hydrogen phosphate	2.5 g
Glucose monohydrate/anhydrous	2.5/2.3 g
Water R	1 000 ml

pH after sterilization 7.1 to 7.5.

68. While sterility testing using membrane filtration, what should be suitable volume of diluent and sample quantity of product to be examined prescribed?

Liquids sample

Minimum quantity to be used for each medium Quantity per container	Minimum quantity to be used for each medium unless otherwise justified and authorized
• less than 1 ml	The whole contents of each container
• 1-40 ml	Half the contents of each container but not less than 1 ml
• greater than 40 ml and not greater than 1	20 ml

X



Insoluble preparations, creams and ointments sample

Minimum quantity to be used for each medium Quantity per container	Minimum quantity to be used for each medium unless otherwise justified and authorized
Insoluble preparations, creams and ointments to be suspended or emulsified	Use the contents of each container to provide not less than 200 mg

Solids sample

Minimum quantity to be used for each medium Quantity per container	Minimum quantity to be used for each medium unless otherwise justified and authorized
• less than 50 mg	The whole contents of each container
• 50 mg or more but less than 300 mg	Half the contents of each container but not less than 50 mg
• 300 mg – 5 g	150 mg
• greater than 5 g	500 mg

Reference: WHO, Document QAS/11.413 FINAL, March 2012



76. What is preferred filter diameter for sterility testing using filtration method?

50 mm in diameter

77. Explain Membrane filtration method for sterility testing?

- The filtration apparatus and membrane are sterilized by appropriate means.
- Use membrane filters with nominal pore size not greater than 0.45 μm and diameter 50 mm.
- Solution to be examined shall be filtered under aseptic conditions.
- Aseptically remove the membrane and transfer to the medium.

Method for aqueous solutions:

- Transfer a small quantity of a suitable, sterile diluent such as a 1 g/l neutral solution of meat or casein peptone pH 6.9 to 7.3 onto the membrane in the apparatus and filter.
- Use diluent containing suitable neutralizing substances and/or inactivating substances in the case of antibiotics.
- Transfer the contents of the container or containers to be tested to the membrane or membranes, if necessary after diluting to the volume used in the method suitability test with the chosen sterile diluent but in any case using not less than the quantities of the product to be examined prescribed in Table-1.



- Do not exceed a washing cycle of five times 100 ml per filter, even if during method suitability it has been demonstrated that such a cycle does not fully eliminate the antimicrobial activity.
- Transfer the whole membrane to the culture medium or cut it aseptically into two equal parts and transfer one half to each of two suitable media.
- Use the same volume of each medium as in the method suitability test.
- Alternatively, transfer the medium onto the membrane in the apparatus. Incubate the media for not less than 14 days.

Soluble solids:

Use for each medium not less than the quantity prescribed in Table 3 of the product dissolved in a suitable solvent such as the solvent provided with the preparation, water for injections R, sodium chloride (9 g/l) TS or peptone (1 g/l) TS1 and proceed with the test as described above for aqueous solutions using a membrane appropriate to the chosen solvent.

Oils and oily solutions:

- Use for each medium not less than the quantity of the product prescribed in Table 2.
- Oils and oily solutions of sufficiently low viscosity may be filtered without dilution through a dry membrane.
- Viscous oils may be diluted as necessary with a suitable sterile diluent such as isopropyl myristate R shown not to have antimicrobial activity in the conditions of the test.
- Allow the oil to penetrate the membrane by its own weight then filter, applying the pressure or suction gradually.
- Wash the membrane at least three times by filtering through it each time about 100 ml of a suitable sterile solution such as peptone (1 g/l) TS1 containing a suitable emulsifying agent at a concentration shown to be appropriate in the method suitability test, for example polysorbate 80 at a concentration of 10 g/l.
- Transfer the membrane or membranes to the culture medium or media or vice versa as described above for aqueous solutions, and incubate at the same temperatures and for the same times.

Ointments and creams:

Table: 3 – Solids sample

Minimum quantity to be used for each medium Quantity per container	Minimum quantity to be used for each medium unless otherwise justified and authorized
• less than 50 mg	The whole contents of each container
• 50 mg or more but less than 300 mg	Half the contents of each container but not less than 50 mg
• 300 mg – 5 g	150 mg
• greater than 5 g	500 mg

Reference: WHO, Document QAS/11.413 FINAL, March 2012

72. Explain the direct inoculation of the culture medium for sterility testing?

- Transfer the quantity of the preparation to be examined prescribed in the following Table directly into the culture medium so that the volume of the product is not more than 10% of the volume of the medium, unless otherwise prescribed.

Table: 1 – Liquids samples

Minimum quantity to be used for each medium Quantity per container	Minimum quantity to be used for each medium unless otherwise justified and authorized
• less than 1 ml	The whole contents of each container
• 1-40 ml	Half the contents of each container but not less than 1 ml
• greater than 40 ml and not greater than 100 ml	20 ml
• greater than 100 ml	10 per cent of the contents of the container but not less than 20 ml
Antibiotic liquids	1 ml

Table: 2 – Insoluble preparations, creams and ointments samples

X

Table: 3 – Solids samples

Minimum quantity to be used for each medium Quantity per container	Minimum quantity to be used for each medium unless otherwise justified and authorized
• less than 50 mg	The whole contents of each container
• 50 mg or more but less than 300 mg	Half the contents of each container but not less than 50 mg
• 300 mg – 5 g	150 mg
• greater than 5 g	500 mg

- If the product to be examined has antimicrobial activity, carry out the test after neutralizing this with a suitable neutralizing substance or by dilution in a sufficient quantity of culture medium.
- When it is necessary to use a large volume of the product it may be preferable to use a concentrated culture medium prepared in such a way that it takes account of the subsequent dilution.

- Where appropriate the concentrated medium may be added directly to the product in its container.

Oily liquids:

- Use media to which have been added a suitable emulsifying agent at a concentration shown to be



/l.

- Transfer the diluted product to a medium not containing an emulsifying agent.
- Incubate the inoculated media for not less than 14 days. Observe the cultures several times during the incubation period.
- Shake cultures containing oily products gently each day. However when fluid thioglycollate medium is used for the detection of anaerobic microorganisms keep shaking or mixing to a minimum in order to maintain anaerobic conditions.**Reference:** WHO, Document QAS/11.413 FINAL, March 2012

[For more interview questions and answers Click Here](#)

[← Previous Post](#)

[Next Post →](#)



[report this ad](#)

Recent Posts

[Courses for Pharmaceutical Industry Professionals](#) January 29, 2023

[Pharma GMP News of the Week January 29, 2023](#) January 29, 2023

[Pharma GMP News of the Week January 22, 2023](#) January 21, 2023

[Pharma GMP News of the Week January 15, 2023](#) January 15, 2023

[Pharma GMP News of the Week January 8, 2023](#) January 7, 2023



Submit your email address to get FREE Weekly Newsletter from <https://tech-publish.com/>

Enter your first name

Enter your last name

Enter your email

Solve This Question

3+3=

Click to Subscribers



Learn anything
Thousands of top courses to choose from

Get started



X


















































































of Amazon.in, or its affiliates.



thinking huts

3D printing access to
education. How we do it.

 ezoic

[report this ad](#)

